

Anti-Cdh19 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H159

Item Number: CM0014929

Introduction

High-quality rabbit polyclonal antibody against human CDH19 (Cadherin-19), validated for WB and ELISA. Suitable for cell adhesion and tissue-sorting studies, with cross-reactivity to human samples.

[Learn More](#)

Parameter	Value
Product Name	CDH19 Rabbit pAb
Remarks/Alias	CDH7; CDH7L2; CDH19
Species	Human
GeneID (Human)	28513
GeneID	28513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-440 of human CDH19 (NP_066976.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H159
Protein Weight	87kDa
Shipping	Ice bag