


Anti-Cdh3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22223

Item Number: CM0017908



Introduction

Rabbit polyclonal CDH3 antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Cadherin-3 (P-cadherin) is a calcium-dependent cell adhesion protein; ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	CDH3 Rabbit pAb
Remarks/Alias	CDHP; HJMD; PCAD; CDH3
Species	Human
GeneID (human)	1001
GeneID	1001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-350 of human CDH3 (NP_001784.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22223
Protein Weight	91kDa
Shipping	Ice bag