

Anti-Cdh4 Rabbit Polyclonal Antibody For Wb, Elisa - P55283

Item Number: CM0016235

Introduction

Rabbit polyclonal antibody targeting CDH4 (cadherin-4, R-cadherin). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Recombinant immunogen, 100 kDa. Suitable for cell adhesion and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	Anti-CDH4 Rabbit Polyclonal Antibody
Remarks/Alias	CAD4; RCAD; R-CAD; CDH4
Species	Human
GeneID (Human)	1002
GeneID	1002
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 757-916 of human CDH4 (NP_001785.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55283
Protein Weight	100kDa
Shipping	Ice bag