

Anti-Cdh8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55286

Item Number: CM0022340

Introduction

Rabbit polyclonal antibody against human Cadherin-8 (CDH8) for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for cell adhesion and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CDH8 Rabbit pAb
Species	Human
GeneID	1006
Immunogen	Recombinant protein of human CDH8
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P55286
Shipping	Ice bag