

Anti-Cdip1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H305

Item Number: CM0016790

Introduction

High-quality rabbit polyclonal antibody against human CDIP1, a p53 apoptosis effector, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	CDIP1 Rabbit pAb
Remarks/Alias	I1; CDIP; LITAF1; C16orf5; CDIP1
Species	Human
GeneID (Human)	29965
GeneID	29965
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CDIP1 (NP_001185983.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H305
Protein Weight	22kDa
Shipping	Ice bag