


Anti-Cdk1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P06493

Item Number: CM0014135



Introduction

High-quality rabbit polyclonal antibody against human CDK1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CDK1 Rabbit pAb
Remarks/Alias	CDC2, CDC28A, P34CDC2, CDK1
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human CDK1 (NP_001777.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag