


Anti-Cdk1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P06493

Item Number: CM0003228



Introduction

CDK1 rabbit monoclonal antibody validated for Western blot, immunoprecipitation, and ELISA. Detects human and mouse CDK1 (P06493), a serine/threonine kinase crucial for cell cycle progression. Suitable for oncology and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CDK1 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; CDK1
Species	Human
GeneID	983
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 198-297 of human CDK1 (NP_001777.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag