


Anti-Cdk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P06493

Item Number: CM0013661



Introduction

High-quality rabbit polyclonal antibody against human CDK1, validated for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat CDK1. Key mitotic kinase for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CDK1 Rabbit pAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; CDK1
Species	Human
Gene ID (Human)	983
Gene ID	983
Immunogen	Synthetic peptide/Recombinant fusion protein containing a sequence corresponding to amino acids 198-297 of human CDK1 (NP_001777.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag