

Anti-Cdk10 Rabbit Polyclonal Antibody For Wb, Elisa - Q15131

Item Number: CM0030336

Introduction

Rabbit polyclonal antibody against human CDK10 (cyclin-dependent kinase 10), validated for WB and ELISA. Detects human and mouse CDK10 (predicted MW 41 kDa). Ideal for cell cycle, cytoskeleton, and ciliogenesis research.

[Learn More](#)

Parameter	Value
Product Name	CDK10 Rabbit pAb
Remarks/Alias	ALSAS; PISSLRE; CDK10
Species	Human
GeneID	8558 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human CDK10 (NP_001092003.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15131
Protein Weight	41kDa
Shipping	Ice bag