

## Anti-Cdk11B Polyclonal Antibody For Wb, If/Icc, Elisa - P21127

Item Number: CM0015214



### Introduction

Rabbit polyclonal antibody against human CDK11B, validated for Western blot, IF/ICC, and ELISA. Reacts with human and mouse. Ideal for studying transcription and splicing regulation, cell cycle, and centromere function.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | CDK11B Rabbit pAb                                                                                                 |
| Remarks/Alias    | p58; PK58; CDK11; CLK-1; CDC2L1; PITSLREA; p58CLK-1; CDK11-p46; CDK11-p58; p58CDC2L1; CDK11-p110; CDK11B          |
| Species          | Human                                                                                                             |
| GeneID           | 984                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human CDK11B (NP_277021.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | P21127                                                                                                            |
| Protein Weight   | 93kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |