

Anti-Cdk11B Rabbit Polyclonal Antibody For If/Icc, Elisa - P21127

Item Number: CM0017860



Introduction

Rabbit polyclonal antibody against CDK11B, validated for IF/ICC and ELISA. Reacts with human, mouse, and rat. Suitable for studying CDK11B's roles in transcription, splicing, and mitosis.

[Learn More](#)

Parameter	Value
Product Name	CDK11B Rabbit pAb
Remarks/Alias	p58; PK58; CDK11; CLK-1; CDC2L1; PITSLREA; p58CLK-1; CDK11-p46; CDK11-p58; p58CDC2L1; CDK11-p110; CDK11B
Species	Human
GeneID (Human)	984
GeneID	984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human CDK11B (NP_277021.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21127
Protein Weight	93kDa
Shipping	Ice bag