


Anti-Cdk12 Monoclonal Antibody For Ihc-P, Elisa - Q9Nyv4

Item Number: CM0006624



Introduction

Anti-CDK12 Rabbit Monoclonal Antibody validated for IHC-P and ELISA applications. Recognizes human, mouse, and rat CDK12, a kinase regulating transcription elongation and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	CDK12 Rabbit mAb
Remarks/Alias	CRK7; CRKR; CRKRS; CDK12
Species	Human
GeneID (Human)	51755
GeneID	51755
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 420-516 of human CDK12 (NP_057591.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYV4
Protein Weight	164kDa
Shipping	Ice bag