


Anti-Cdk12 Rabbit Pab For Wb, Ihc-P, Elisa - Q9Nyv4

Item Number: CM0018353



Introduction

Rabbit polyclonal CDK12 antibody targeting the N-terminal region (aa 1-300). Applications: WB, IHC-P, ELISA. Reactivity: human, mouse, rat. Ideal for studying transcription elongation and genomic stability.

[Learn More](#)

Parameter	Value
Product Name	CDK12 Rabbit pAb
Remarks/Alias	CRK7, CRKR, CRKRS, CDK12
Species	Human
Gene ID (Human)	51755
Gene ID	51755
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CDK12 (NP_057591.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYV4
Protein Weight	164kDa
Shipping	Ice bag