

Anti-Cdk16 Polyclonal Antibody For If/Icc, Elisa - Q00536

Item Number: CM0024831



Introduction

Rabbit polyclonal antibody against human CDK16 (PCTAIRE1, PCTK1) for IF/ICC and ELISA. Validated for cross-reactivity with human, target weight 56kDa.

[Learn More](#)

Parameter	Value
Product Name	CDK16 Rabbit pAb
Remarks/Alias	PCTK1; PCTAIRE; PCTAIRE1; PCTGAIRE; CDK16
Species	Human
Gene ID (Human)	5127
Gene ID	5127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 277-496 of human CDK16 (NP_006192.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q00536
Protein Weight	56kDa
Shipping	Ice bag