

Anti-Cdk18 Polyclonal Antibody For Wb And Elisa - Q07002

Item Number: CM0023398

Introduction

Rabbit polyclonal antibody targeting human CDK18 (PCTK3/PCTAIRE3). Recombinant immunogen (aa 373-474). Validated for Western blot and ELISA. Cross-reacts with human and mouse. Ideal for signal transduction studies in differentiated cells.

[Learn More](#)

Parameter	Value
Product Name	CDK18 Rabbit pAb
Remarks/Alias	PCTK3, PCTAIRE, PCTAIRE3, CDK18
Species	Human
GeneID (human)	5129
GeneID	5129
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 373-474 of human CDK18 (NP_002587.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q07002
Protein Weight	54kDa
Shipping	Ice bag