

Anti-Cdk19 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bwu1

Item Number: CM0023741

Introduction

Rabbit polyclonal antibody targeting human CDK19, a kinase component of the Mediator complex involved in transcriptional regulation. Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for transcription and kinase research.

[Learn More](#)

Parameter	Value
Product Name	CDK19 Rabbit pAb
Remarks/Alias	CDK11; DEE87; CDC2L6; EIEE87; bA346C16.3; CDK19
Species	Human
Gene ID (Human)	23097
Gene ID	23097
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CDK19 (XP_005266928.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BWU1
Protein Weight	57kDa
Shipping	Ice bag