


Anti-Cdk2 (Phospho-T160) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P24941

Item Number: CM0001781



Introduction

Phospho-CDK2-T160 Rabbit mAb validated for WB, IF-P, IHC-P, ELISA; reacts with human, mouse, rat CDK2; targets phosphorylated T160; ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK2-T160 Rabbit mAb
Remarks/Alias	CDKN2; p33(CDK2); Phospho-CDK2-T160
Species	Human
GeneID (Human)	1017
GeneID	1017
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T160 of human CDK2 (NP_001789.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24941
Protein Weight	34kDa
Shipping	Ice bag