


Anti-Cdk2 Rabbit Monoclonal Antibody (Ko Validated) For Western Blot And Elisa - P24941

Item Number: CM0003252



Introduction

Anti-CDK2 rabbit monoclonal antibody (KO validated) for Western blot and ELISA. Confirmed reactivity with human, mouse, and rat CDK2. Validated for knockout studies. Targets the cyclin-dependent kinase 2 (P24941), a cell cycle master regulator.

[Learn More](#)

Parameter	Details
Product Name	[KO Validated] CDK2 Rabbit mAb
Remarks/Alias	CDKN2; p33(CDK2); K2
Species	Human
GeneID (Human)	1017
GeneID	1017
Immunogen	Synthetic Peptide corresponding to amino acids 199-298 of human CDK2 (P24941)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24941
Protein Weight	30 kDa/34 kDa
Shipping	Ice bag