


Anti-Cdk2 Rabbit Polyclonal Antibody [Ko Validated] For Wb, Ihc-P, Elisa - P24941

Item Number: CM0022176



Introduction

Rabbit polyclonal antibody against human CDK2, validated for WB, IHC-P, and ELISA. Detects endogenous CDK2 (~30/34 kDa) in human, mouse, and rat samples. Ideal for cell cycle, DNA repair, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK2 Rabbit pAb
Remarks/Alias	CDKN2; p33(CDK2); K2
Species	Human
GeneID (Human)	1017
GeneID	1017
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-298 of human CDK2 (NP_001789.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24941
Protein Weight	30 kDa/34 kDa
Shipping	Ice bag