

Anti-Cdk20 Polyclonal Antibody For Ihc-P And Elisa - Q8Izl9

Item Number: CM0025559

Introduction

CDK20 rabbit polyclonal antibody validated for IHC-P and ELISA, reacting with human and rat. Ideal for cell cycle, Shh signaling, and cilogenesis research applications.

[Learn More](#)

Parameter	Value
Product Name	CDK20 Rabbit pAb
Remarks/Alias	P42; CCRK; CDCH; PNQALRE; CDK20
Species	Human
GeneID (Human)	23552
GeneID	23552
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human CDK20 (Q8IZL9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Rat
SWISS	Q8IZL9
Protein Weight	39kDa
Shipping	Ice bag