


Anti-Cdk2Ap1 Rabbit Polyclonal Antibody For Wb, Elisa - O14519

Item Number: CM0020379



Introduction

CDK2AP1 Rabbit Polyclonal Antibody (pAb) validated for Western blot and ELISA. Targets human CDK2AP1 (DOC-1), a CDK2 inhibitor and component of the NuRD chromatin remodeling complex. Ideal for cancer and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	CDK2AP1 Rabbit Polyclonal Antibody
Remarks/Alias	DOC1; ST19; DORC1; doc-1; p12DOC-1; CDK2AP1
Species	Human
Gene ID (Human)	8099
Gene ID	8099
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human CDK2AP1 (NP_004633.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14519
Protein Weight	12kDa
Shipping	Ice bag