

Anti-Cdk4 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P11802

Item Number: CM0017007



Introduction

CDK4 rabbit polyclonal antibody validated for WB, IF/ICC, IP, and ELISA. Reacts with human and mouse CDK4, a key cell cycle kinase regulating G1/S transition. Ideal for cancer research and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CDK4 Rabbit pAb
Remarks/Alias	CMM3; PSK-J3; CDK4
Species	Human
Gene ID (Human)	1019
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-303 of human CDK4 (NP_000066.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11802
Protein Weight	34kDa
Shipping	Ice bag