

Anti-Cdk4 Rabbit Polymab® For Wb, Ip, If/Icc, If-P, Ihc-P, Elisa - P11802

Item Number: CM0004311



Introduction

Rabbit monoclonal antibody targeting human CDK4, suitable for WB, IP, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK4 Rabbit PolymAb®
Remarks/Alias	CMM3; PSK-J3
Species	Human
GeneID (human)	1019
GeneID	1019
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 204-303 of human CDK4 (NP_000066.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11802
Protein Weight	34kDa
Shipping	Ice bag