


Anti-Cdk5 Rabbit Monoclonal Antibody [Ko Validated] For Wb And Elisa - Q00535

Item Number: CM0006925



Introduction

KO-validated rabbit monoclonal antibody targeting CDK5, suitable for WB and ELISA. Reacts with human, mouse, and rat. Explores neuronal development, apoptosis, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK5 Rabbit mAb
Remarks/Alias	LIS7; PSSALRE; K5
Species	Human
Gene ID (Human)	1020
Gene ID	1020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 92-289 of human CDK5 (NP_004926.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00535
Protein Weight	33kDa
Shipping	Ice bag