


Anti-Cdk5 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q00535

Item Number: CM0028533



Introduction

Validated CDK5 rabbit polyclonal antibody for WB, IF/ICC, IP, ELISA. Recognizes human, mouse, rat CDK5. Essential for neuronal development. Ideal for research on cell cycle, apoptosis, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK5 Rabbit pAb
Remarks/Alias	LIS7; PSSALRE; K5
Species	Human
GeneID (Human)	1020
GeneID	1020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human CDK5 (NP_004926.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00535
Protein Weight	33kDa
Shipping	Ice bag