


Anti-Cdk5 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - Q00535

Item Number: CM0022581



Introduction

KO validated rabbit polyclonal antibody against human CDK5 (cyclin-dependent kinase 5), suitable for Western blot (WB) and ELISA. Detects CDK5 in human, mouse, and rat samples, supporting research in cell cycle, neurogenesis, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK5 Rabbit pAb
Remarks/Alias	LIS7; PSSALRE; K5
Species	Human
Gene ID	1020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-289 of human CDK5 (NP_004926.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00535
Protein Weight	33kDa
Shipping	Ice bag