

Anti-Cdk5R1 Monoclonal Antibody For Wb, Elisa - Q15078

Item Number: CM0006809



Introduction

CDK5R1 Rabbit mAb (p35) validated for human, mouse, rat WB/ELISA. Neuron-specific CDK5 activator involved in neurite outgrowth, cortical lamination, and circadian clock regulation.

[Learn More](#)

Parameter	Value
Product Name	CDK5R1 Rabbit mAb
Remarks/Alias	p23; p25; p35; CDK5R; NCK5A; CDK5P35; p35nck5a; CDK5R1
Species	Human
GeneID	8851
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human CDK5R1 (NP_003876.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15078
Protein Weight	34kDa
Shipping	Ice bag