


Anti-Cdk5R1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q15078

Item Number: CM0007597



Introduction

High-quality rabbit monoclonal antibody against CDK5R1 (p35) suitable for Western Blot, IF/ICC, IF-P, and ELISA in Human, Mouse, and Rat. Ideal for neuroscience research on neurite outgrowth and circadian rhythms.

[Learn More](#)

Parameter	Value
Product Name	CDK5R1 Rabbit mAb
Remarks/Alias	p23; p25; p35; CDK5R; NCK5A; CDK5P35; p35nck5a; CDK5R1
Species	Human
GeneID (Human)	8851
GeneID	8851
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15078
Protein Weight	34kDa
Shipping	Ice bag