


## Anti-Cdk5Rap2 Polyclonal Antibody For Wb, Elisa - Q96Sn8

Item Number: CM0016405



### Introduction

Rabbit polyclonal antibody targeting CDK5RAP2 (Cep215) for Western blot and ELISA applications. Reacts with human; 215 kDa. For studies on centriole regulation and spindle orientation.

[Learn More](#)

| Specification    | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CDK5RAP2 Rabbit pAb                                                                                                      |
| Remarks/Alias    | C48; MCPH3; Cep215; CDK5RAP2                                                                                             |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 55755                                                                                                                    |
| GeneID           | 55755                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1744-1893 of human CDK5RAP2 (NP_060719.4). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, ELISA                                                                                                                |
| Cross Reactivity | Human                                                                                                                    |
| SWISS            | Q96SN8                                                                                                                   |
| Protein Weight   | 215kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |