


Anti-Cdk5Rap3 (C53/Lzap) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Jb5

Item Number: CM0027852



Introduction

Rabbit polyclonal antibody against CDK5RAP3 (C53/LZAP), raised with recombinant immunogen (aa 1-300). Validated in WB, IF/ICC, ELISA for human, mouse, rat. Suitable for studies on ufmylation, ribosome recycling, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	CDK5RAP3 Rabbit pAb
Remarks/Alias	C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114; CDK5RAP3
Species	Human
GeneID (Human)	80279
GeneID	80279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CDK5RAP3 (NP_788276.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96JB5
Protein Weight	57kDa
Shipping	Ice bag