


Anti-Cdk6 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - Q00534

Item Number: CM0009042



Introduction

Rabbit monoclonal antibody specific to human CDK6 (37 kDa), a key regulator of G1/S cell cycle transition. Validated for WB, IP, IF/ICC, and ELISA. Useful for oncology and cell differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	CDK6 Rabbit mAb
Remarks/Alias	MCPH12; PLSTIRE
Species	Human
GeneID (human)	1021
GeneID	1021
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 227-326 of human CDK6 (NP_001250.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q00534
Protein Weight	37 kDa
Shipping	Ice bag