


Anti-Cdk6 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q00534

Item Number: CM0013692



Introduction

Rabbit polyclonal antibody against human CDK6, suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse. Detects a band at ~37 kDa. Ideal for cell cycle regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CDK6 Rabbit pAb
Remarks/Alias	MCPH12; PLSTIRE; CDK6
Species	Human
Gene ID (Human)	1021
Gene ID	1021
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 226-326 of human CDK6 (NP_001250.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00534
Protein Weight	37kDa
Shipping	Ice bag