


Anti-Cdk6 Rabbit Polyclonal Antibody For Wb - Q00534

Item Number: CM0014497



Introduction

CDK6 Rabbit pAb validated for WB in human, mouse, rat. Targets cyclin-dependent kinase 6 (37kDa), a key cell cycle regulator. Ideal for cancer and differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	CDK6 Rabbit pAb
Remarks/Alias	MCPH12; PLSTIRE; CDK6
Species	Human
GeneID	1021 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-326 of human CDK6 (NP_001250.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00534
Protein Weight	37kDa
Shipping	Ice bag