


Anti-Cdk8 Rabbit Monoclonal Antibody For Wb And Elisa - P49336

Item Number: CM0011616



Introduction

Rabbit monoclonal antibody targeting human cyclin-dependent kinase 8 (CDK8). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for Mediator complex and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CDK8 Rabbit mAb
Remarks/Alias	K35; IDDHBA; CDK8
Species	Human
GeneID (Human)	1024
GeneID	1024
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human CDK8 (P49336).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P49336
Protein Weight	53kDa
Shipping	Ice bag