


Anti-Cdk8 Rabbit Monoclonal Antibody For Wb And Elisa - P49336

Item Number: CM0005382



Introduction

High specificity anti-CDK8 rabbit monoclonal antibody validated for Western blot and ELISA. Recognizes human and mouse CDK8, a kinase essential for Mediator complex function and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	CDK8 Rabbit mAb
Remarks/Alias	CDK8; K35; cyclin-dependent kinase 8
Species	Human
GeneID	1024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 364-464 of human CDK8 (NP_001251.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49336
Protein Weight	53kDa
Shipping	Ice bag