

Anti-Cdk9 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Ip, Elisa - P50750

Item Number: CM0003075



Introduction

KO validated rabbit monoclonal anti-CDK9 antibody detects human, mouse, rat CDK9. Suitable for WB, IHC-P, IF/ICC, IP, ELISA. Ideal for transcription and gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK9 Rabbit mAb
Remarks/Alias	TAK; C-2k; CTK1; CDC2L4; PITALRE; K9
Species	Human
GeneID (Human)	1025
GeneID	1025
Immunogen	Synthetic peptide corresponding to amino acids 273-372 of human CDK9 (P50750)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50750
Protein Weight	43kDa
Shipping	Ice bag