


Anti-Cdk9 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50750

Item Number: CM0013666



Introduction

Human CDK9 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat CDK9. Suitable for transcription regulation and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	CDK9 Rabbit pAb
Remarks/Alias	TAK; C-2k; CTK1; CDC2L4; PITALRE; CDK9
Species	Human
GeneID (human)	1025
GeneID	1025
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50750
Protein Weight	43kDa
Shipping	Ice bag