

Anti-Cdkl3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Ivw4

Item Number: CM0019603



Introduction

Rabbit polyclonal antibody against CDKL3 (NKIAMRE) validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for kinase signaling research. UniProt Q8IVW4.

[Learn More](#)

Parameter	Value
Product Name	CDKL3 Rabbit pAb
Remarks/Alias	NKIAMRE; CDKL3
Species	Human
Gene ID (Human)	51265
Gene ID	51265
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human CDKL3 (NP_001107047)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IVW4
Protein Weight	68kDa
Shipping	Ice bag