

Anti-Cdkn1A Rabbit Polyclonal Antibody For Wb, Elisa - P39689

Item Number: CM0022230

Introduction

High-quality anti-CDKN1A rabbit polyclonal antibody validated for WB and ELISA. Reacts with mouse and rat CDKN1A (p21), a key cell cycle inhibitor. Ideal for cancer and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	CDKN1A Rabbit pAb
Remarks/Alias	P21; CDKI; CIP1; SDI1; Waf1; mda6; CAP20; Cdkn1; p21WAF; p21Cip1; CDKN1A
Species	Mouse
GeneID (human)	1026
GeneID	12575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-159 of mouse CDKN1A (NP_031695.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P39689
Protein Weight	18kDa
Shipping	Ice bag