

Anti-Cdkn1A/P21Cip1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P38936

Item Number: CM0007224



Introduction

Rabbit monoclonal antibody against human CDKN1A/p21CIP1 (p21). Validated for Western blot, immunohistochemistry, and ELISA. Detects an 18 kDa band; ideal for cell cycle and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	CDKN1A/p21CIP1 Rabbit mAb
Remarks/Alias	P21; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1; CDKN1A/p21CIP1
Species	Human
GeneID (Human)	1026
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human CDKN1A/p21CIP1 (NP_000380.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P38936
Protein Weight	18kDa
Shipping	Ice bag