

Anti-Cdkn1A/P21Cip1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - P38936

Item Number: CM0012728



Introduction

Rabbit monoclonal antibody against human CDKN1A/p21CIP1, validated for WB, IP, IF/ICC, ELISA. Recognizes human and mouse targets; ideal for cell cycle and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDKN1A/p21CIP1 Rabbit mAb
Remarks/Alias	P21; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1
Species	Human
GeneID (Human)	1026
GeneID	1026
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human CDKN1A/p21CIP1 (NP_000380.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P38936
Protein Weight	18kDa
Shipping	Ice bag