

Anti-Cdkn1A/P21Cip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P38936

Item Number: CM0015941



Introduction

High-quality rabbit polyclonal antibody against CDKN1A/p21CIP1 for Western blot, IHC, IF, and ELISA. Validated in human, mouse, rat. Suitable for cell cycle and DNA damage response research.

[Learn More](#)

Parameter	Value
Product Name	CDKN1A/p21CIP1 Rabbit pAb
Remarks/Alias	P21; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1; CDKN1A/p21CIP1
Species	Human
GeneID (Human)	1026
GeneID	1026
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-164 of human CDKN1A/p21CIP1 (NP_000380.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38936
Protein Weight	18kDa
Shipping	Ice bag