


Anti-Cdkn1B/P27Kip1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P46527

Item Number: CM0017717



Introduction

Rabbit polyclonal antibody against CDKN1B/p27KIP1, validated for WB, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat. Cyclin-dependent kinase inhibitor regulating cell cycle progression; essential for G1 arrest. Suitable for tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	CDKN1B/p27KIP1 Rabbit pAb
Remarks/Alias	KIP1; MEN4; CDKN4; MEN1B; P27KIP1; CDKN1B/p27KIP1
Species	Human
GeneID (Human)	1027
GeneID	1027
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human CDKN1B/p27KIP1 (NP_004055.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46527
Protein Weight	22kDa
Shipping	Ice bag