


Anti-Cdkn1B/P27Kip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P46527

Item Number: CM0020067



Introduction

High-quality rabbit polyclonal antibody targeting human CDKN1B/p27KIP1 (P46527). Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Suitable for cell cycle regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CDKN1B/p27KIP1 Rabbit pAb
Remarks/Alias	KIP1; MEN4; CDKN4; MEN1B; P27KIP1; CDKN1B/p27KIP1
Species	Human
Gene ID (Human)	1027
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-198 of human CDKN1B/p27KIP1 (NP_004055.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46527
Protein Weight	22kDa
Shipping	Ice bag