

Anti-Cdkn1C/P57 Kip2 Polyclonal Antibody For Wb, Elisa - P49918

Item Number: CM0023382

Introduction

Rabbit polyclonal antibody against CDKN1C/p57 Kip2 (cyclin-dependent kinase inhibitor 1C), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	CDKN1C/p57 Kip2 Rabbit pAb
Remarks/Alias	BWS; WBS; p57; BWCR; KIP2; p57Kip2; CDKN1C/p57 Kip2
Species	Human
Gene ID (Human)	1028
Gene ID	1028
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human CDKN1C/p57 Kip2 (NP_000067.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49918
Protein Weight	32kDa
Shipping	Ice bag