

Anti-Cdkn1C/P57 Kip2 Rabbit Monoclonal Antibody For Wb - P49918

Item Number: CM0005958

Introduction

Rabbit monoclonal antibody against human CDKN1C/p57 Kip2 (cyclin-dependent kinase inhibitor 1C, P49918), validated for Western blotting and ELISA. Ideal for studying cell cycle arrest, proliferation, and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	CDKN1C/p57 Kip2 Rabbit mAb
Remarks/Alias	BWS; WBS; p57; BWCR; KIP2; p57Kip2; CDKN1C/p57 Kip2
Species	Human
Gene ID (Human)	1028
Gene ID	1028
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49918
Protein Weight	32kDa
Shipping	Ice bag