


Anti-Cdkn2A/P16Ink4A Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P42771

Item Number: CM0003255



Introduction

High-quality [KO Validated] rabbit monoclonal antibody targeting human CDKN2A/p16INK4a (SWISS P42771). Validated for WB, IF/ICC, and ELISA, suitable for cell cycle and tumor suppressor research applications.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDKN2A/p16INK4a Rabbit mAb
Remarks/Alias	ARF, MLM, P14, P16, P19, CMM2, INK4, MTS1, TP16, CDK4I, CDKN2, INK4A, MTS-1, P14ARF, P19ARF, P16INK4, P16INK4A, P16-INK4A, 4a
Species	Human
GeneID	1029
Immunogen	Synthetic Peptide immunogen. A synthetic peptide corresponding to a sequence within amino acids 57-156 of human CDKN2A/p16INK4a (P42771).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P42771
Protein Weight	17 kDa
Shipping	Ice bag