


Anti-Cdkn2A/P19Arf Rabbit Monoclonal Antibody For Wb, Elisa - Q64364

Item Number: CM0007493



Introduction

CamelBio offers anti-mouse CDKN2A/p19ARF rabbit monoclonal antibody validated for WB and ELISA. Targets tumor suppressor ARF, suitable for cell cycle, apoptosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CDKN2A/p19ARF Rabbit mAb
Remarks/Alias	Arf; p16; MTS1; Pctr1; p19ARF; p16INK4a; p19<ARF>; ARF-INK4a; INK4a-ARF; Ink4a/Arf; p16(INK4a); CDKN2A/p19ARF
Species	Mouse
GeneID (human)	1029
GeneID	12578
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q64364
Protein Weight	19KDa
Shipping	Ice bag