

Anti-Cdkn2A/P19Arf Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q64364-1

Item Number: CM0006188

Introduction

Rabbit monoclonal antibody against mouse CDKN2A/p19ARF (p19ARF). For WB, IF/ICC, IF-P, ELISA. Immunogen: synthetic peptide within amino acids 1-100. Cross-reacts with mouse. Calculated molecular weight: 19 kDa.

[Learn More](#)

Parameter	Value
Product Name	CDKN2A/p19ARF Rabbit mAb
Remarks/Alias	Arf; p16; MTS1; Pctr1; p19ARF; p16INK4a; p19<ARF>; ARF-INK4a; INK4a-ARF; Ink4a/Arf; p16(INK4a); CDKN2A/p19ARF
Species	Mouse
Gene ID (Human)	1029
Gene ID	12578
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse p19ARF (NP_034007.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q64364-1
Protein Weight	19 kDa
Shipping	Ice bag