

Anti-Cdkn2B Rabbit Polyclonal Antibody For Wb, Elisa - P42772

Item Number: CM0026581

Introduction

CDKN2B Rabbit Polyclonal Antibody (P42772) validated for WB and ELISA, reactive with human, mouse, and rat. Ideal for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Specification
Product Name	CDKN2B Rabbit pAb
Remarks/Alias	P15; MTS2; TP15; CDK4I; INK4B; p15INK4b; CDKN2B
Species	Human
GeneID (Human)	1030
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-138 of human CDKN2B (NP_004927.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42772
Protein Weight	15kDa
Shipping	Ice bag